EN**NASO LIBERO®****Hypertonic Nasal Spray with Aloe and Chamomile****RECOMENDED USE**

For daily nasal hygiene.

A nasal irrigation spray to help relieve nasal congestion associated with sinusitis, blocked nose, hay fever, allergies and post-operative nasal care. Can be used by adults and children from 3 years of age.

NASO LIBERO Hypertonic nasal spray is a non-medicated, preservative free, deeply penetrating nasal irrigation spray and can be used by children and adults. It works by thinning mucus and washing away foreign particles such as dust, pollen and allergens. This helps relieve congestion and restore normal respiratory function. NASO LIBERO Hypertonic nasal spray is non habit forming and can be used together with other medications. Frequent sufferers of rhinitis and sinusitis may find that regular use of NASO LIBERO Hypertonic nasal spray may help decrease the use of medicated products.

The addition of Aloe and Chamomile helps soothing sore noses through hydrating and soothing properties as well as adding a pleasant fresh aromatic touch.

DIRECTIONS

1. Remove the over cap.
2. Keep head upright and spray desired dose into each nostril. Repeat as often as necessary.

Dosage depends on how long the actuator is activated.

COMPOSITION

Buffered Hypertonic Sea salt solution (2.1%), with addition of 0.025% Aloe Vera and 0.025% Chamomile

WARNING

- The use of this device by more than one person may spread infections, clean nozzle in soapy water and rinse with hot water after each use
- Because of its content, this spray shall not be used if you are hypersensitive (allergic) to chamomile or aloe extract
- Pressurized container: May burst if heated
- Keep away from heat / sparks / open flames / hot surfaces. — No smoking
- Pressurized container: Do not pierce or burn, even after use
- Protect from sunlight
- Do not expose to temperatures exceeding 50°C (122°F)
- Keep out of reach of children

125 ml

Look for other **NASO LIBERO®** products

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Ver: 2015-10-28

NASO LIBERO Nasal spray to be used from 6 months

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Aurena Laboratories AB